



Flow cytometry

Let's get to the science



Analyze single cells like never before

Within a single cell are insights waiting to be translated, answers waiting to be found, and discoveries waiting to be made. Get to them in new ways, in this new era of flow cytometry.

Harness the power of flow at thermofisher.com/flowcytometry



For Research Use Only. Not for use in diagnostic procedures. © 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. COL22153 0917

Simplifying Flow Cytometry

Thermo Fisher Scientific Inc. is the world leader in serving science, with revenues of more than \$20 billion and approximately 65,000 employees globally. Our mission is to enable our customers to make the world healthier, cleaner and safer. We help our customers accelerate life sciences research, solve complex analytical challenges, improve patient diagnostics, deliver medicines to market and increase laboratory productivity.

Through our premier brands – Thermo Scientific, Applied Biosystems, Invitrogen, Fisher Scientific and Unity Lab Services – we offer an unmatched combination of innovative technologies, purchasing convenience and comprehensive services.

For more information, please visit www.thermofisher.com

Harness the power of flow at thermofisher.com/flowcytometry



No-Wash No-Lyse Techniques: Simplifying Flow Cytometry

Conventional flow cytometry has relatively limited sensitivity, meaning that heterogeneous cell samples often have to be purified through lysis and wash steps prior to analysis. Technological advancements have enabled “no-wash, no-lyse” workflows for more precise and faster flow cytometry analysis.



No-Wash/No-Lyse Techniques

Sponsored by: **ThermoFisher**
SCIENTIFIC

Custom publishing from: **TheScientist**
EXPLORING LIFE, INSPIRING INNOVATION



Bound to work



We're testing Invitrogen antibodies to meet the newest specificity standards

We're committed to superior antibody performance. That's why we're rigorously testing our Invitrogen™ antibodies using a two-part approach for advanced verification—so you can be confident your Invitrogen antibody will bind to the correct target in your intended application.

Learn about our two-part approach to antibody testing at thermofisher.com/antibodyvalidation



For Research Use Only. Not for use in diagnostic procedures. © 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. The use or any variation of the word “validation” refers only to research use antibodies that were subject to functional testing to confirm that the antibody can be used with the research techniques indicated. It does not ensure that the product(s) was validated for clinical or diagnostic uses. COL22153 0917

No-Wash/No-Lyse Techniques: Simplifying Flow Cytometry

Conventional flow cytometry has relatively limited sensitivity, meaning that heterogeneous cell samples often have to be purified through lysis and wash steps prior to analysis. Technological advancements have enabled no-wash/no-lyse workflows for more precise and faster flow cytometry analysis.

What are No-Wash/No-Lyse Techniques?

Wash and lysis steps remove unwanted cells from heterogeneous cell samples, leaving purified populations for analysis. However, these steps can cause significant collateral damage to investigated cells [1-3], which contributes to interexperimental and interlaboratory variation [4].

Improved technology has resulted in no-wash/no-lyse protocols for detecting cells without labor-intensive and potentially deleterious wash- or lysis-mediated purification. No-wash/no-lyse workflows have been used to generate accurate, reproducible results for a variety of applications [2,4-6].



No-Wash/No-Lyse Applications


No-wash/no-lyse protocols have become popular for a range of applications in a variety of fields.

Used for:

- Determining absolute cell counts
- Human blood diagnostic analysis [4,6]
- Immunophenotyping [11]
- Identifying changes in protein expression/phosphorylation [2] or enzyme activity/interactions [5]

Used in:

- Vaccine development
- Longitudinal studies
- Oncology
- Cardiovascular disease



NO-WASH/NO-LYSE WORKFLOW



TRADITIONAL WORKFLOW

A Simpler and Faster Workflow

No-wash/no-lyse protocols present researchers with simpler, faster, and more accurate options. Removing wash and lysis steps can reduce workflow duration significantly. Fewer manipulations help to preserve cell integrity, helping promote data integrity.


The Importance of Validated Antibodies

Antibody-based tests and protocols are prevalent at all levels of research. To generate trustworthy and accurate data, antibodies need to be specific, selective, and reproducible. However, no universally accepted best practice guidelines for antibody validation currently exist [9,10].

Antibody non-specificity and non-reproducibility pose threats to data integrity and waste time, energy, and resources. Proper antibody validation for both specificity and functionality is critical.

Vendor-supplied validation information can range broadly [9]. The best vendors will:

- Provide detailed descriptions of validation procedures
- Use proper controls
- Test multiple applications and conditions
- Test multiple antigen sources



Technology Marches On

No-wash/No-lyse workflows take advantage of new flow cytometry technologies and techniques.

Acoustic Cytometry

Acoustic radiation pressure can focus particles into a single tight line for flow cytometric optical interrogation. This acoustic focusing eliminates the need for sample dilution and provides superior flow control compared to conventional hydrodynamic focusing, resulting in greater sensitivity, precision, throughput, and applicative range [7,8].

True Volumetric Absolute Counting (TVAC)

TVAC is based on the principle that $c = N / V$, where the particle concentration (c) is equal to the counted number (N) divided by the volume (V). TVAC directly measures V, yielding an absolute value and providing better precision than conventional volume assessment using calibration counting beads [1].

REFERENCES

1. B. Grove et al., "A New No-Lyse, No-Wash Flow-Cytometric Method for the Determination of CD4+ T Cells in Blood Samples," *Transfus Med Hemother* 30: 8-13, 2003.
2. J. Demaret et al., "Intracellular Flow Cytometry Improvements in Clinical Studies," *Methods Mol Biol* 1524: 315-327, 2017.
3. B. Grove et al., "High-Grade Loss of Leukocytes and Hematopoietic Progenitor Cells Caused by Erythrocyte-Lysing Procedures for Flow Cytometric Analyses," *J Hematother Stem Cell Res* 12(3): 321-330, 2003.
4. A. Alvarez-Larran et al., "A multicolor, no-lyse no-wash assay for the absolute counting of CD34+ cells by flow cytometry," *Cytometry* 50: 247-253, 2002.
5. S. Balaguer et al., "Real-Time Cytometric Assay of Nitric Oxide and Superoxide Interaction in Peripheral Blood Monocytes: A No-Wash, No-Lyse Kinetic Method," *Cytometry Part B* 92B: 211-217, 2017.
6. YX. Cui et al., "Stepwise optimization of the procedure for assessment of circulating progenitor cells in patients with myocardial infarction," *PLoS One* 7(1): e30399, 2012.
7. M.E. Pfyfer et al., "Multicolor acoustic focusing for parallel flow cytometry," *Anal Chem* 84(4): 1831-1839, 2012.
8. M. Ward et al., "Fundamentals of Acoustic Cytometry," *Curr Protoc Cytom* 52:1.22.1-1.22.12, 2009.
9. J. Bordeaux et al., "Antibody validation," *BioTechniques* 48(3): 197-209, 2010.
10. M. Uhlen et al., "A proposal for validation of antibodies," *Acta Med Scand* 131(10): 823-827, 2016.
11. H. Jandis et al., "Lymphocyte and monocyte flow cytometry immunophenotyping as a diagnostic tool in uncharacteristic inflammatory disorders," *BMC Infect Dis* 10:205, 2010.